

Analytical Performance of an Inductively Coupled Plasma Optical Emission Spectrometry with Dual View Configuration

Júlio C. J. Silva, Nivaldo Baccan and Joaquim A. Nóbrega
2003, vol. 14, No 2, 210-315.

In the Table 2 of the page 311, the correct values of IP and EP to Ba II (455.4 nm), Al I (308.2 nm), Ba II (233.5 nm), Mn II (260.6 nm), Mn II (257.6 nm), Mo II (203.8 nm), Mo II (202.0 nm) are 5.21 and 2.72 eV, 5.99 and 4.02 eV,

5.21 and 6.0, 7.44 and 4.75, 7.44 and 4.81, 7.1 and 6.08, 7.1 and 6.13, respectively.

In the Table 2 of page 311, the correct value of EP to Mg I (285.2 nm) is 4.34 eV.